INSTRUCTIONS FOR USE

INTENDED USE

Cell-Free DNA BCT® CE is a direct draw whole blood collection tube intended for collection, stabilization and transportation of cell-free plasma DNA. This device also stabilizes and preserves cellular genomic DNA present in nucleated blood cells and circulating epithelial cells (tumor cells) found in whole blood. **This product is FOR EXPORT ONLY, not to be sold in the United States.**

SUMMARY AND PRINCIPLES

Accurate analysis of cf-DNA can be compromised by sample handling, shipping and processing, causing lysis of nucleated blood cells and subsequent release of cellular genomic DNA. Additionally, degradation of cf-DNA due to blood activity can be problematic.

The formaldehyde-free preservative reagent contained in Cell-Free DNA BCT CE^{1,2} stabilizes nucleated blood cells, preventing the release of cellular genomic DNA, and inhibits nuclease mediated degradation of cf-DNA, contributing to the overall stabilization of cf-DNA³. Samples collected in Cell-Free DNA BCT CE are stable for up to 14 days at temperatures between 6-37° C, allowing convenient sample collection, transport and storage⁴.

The formaldehyde-free preservative reagent contained in Cell-Free DNA BCT CE stabilizes circulating epithelial cells (tumor cells) in whole blood for up to 4 days at temperatures between 15-30° C⁵.

REAGENTS

Cell-Free DNA BCT CE contains the anticoagulant K₂EDTA and a cell preservative in a liquid medium.

PRECAUTIONS

- For In Vitro Diagnostic Use.
- Do not freeze specimens collected in Cell-Free DNA BCT CE as breakage could result.
- Do not use tubes after expiration date.
- Do not use tubes for collection of materials to be injected into patients.
- Product is intended for use as supplied. Do not dilute or add other components to Cell-Free DNA BCT CE.
- Overfilling or underfilling of tubes will result in an incorrect blood-to-additive ratio and may lead to incorrect analytic results or poor product performance.
- CAUTION**
 - Glass has the potential for breakage; precautionary measures should be taken during handling.
 - All biological specimens and materials coming in contact with them are considered biohazards and should be treated as if capable of transmitting infection. Dispose of in accordance with federal, state and local regulations. Avoid contact with skin and mucous membranes.
 - Product should be disposed with infectious medical waste.
 - Remove stopper by either gently rocking the stopper from side to side or by grasping with a simultaneous twisting and pulling action. A "thumb roll" procedure for stopper removal is not recommended, as tube breakage and injury may result. Reinsert stopper by gently pushing stopper onto tube with a simultaneous twisting action.
- SDS can be obtained at www.streck.com or by calling 402-333-1982.

STORAGE AND STABILITY


- When stored at 18-30° C, unused Cell-Free DNA BCT CE is stable through expiration date.
- Do not freeze unfilled Cell-Free DNA BCT CE. Proper insulation may be required for shipment during extreme temperature conditions.
- Blood samples collected in Cell-Free DNA BCT CE for cf-DNA analysis are stable for 14 days when stored between 6-37° C.
- Blood samples collected in Cell-Free DNA BCT CE for genomic DNA analysis are stable for 14 days when stored between 6-37° C.
- Blood samples collected in Cell-Free DNA BCT CE for circulating epithelial cells (tumor cells) are stable for 4 days when stored between 15°-30° C.

INDICATIONS OF PRODUCT DETERIORATION

- Cloudiness or precipitate visible in reagent of unused tube.
- If indications of product deterioration occur, contact Streck Technical Services at 402-691-7510 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

- Collect specimen by venipuncture according to CLSI H3-A6⁶.
Prevention of Backflow - Since Cell-Free DNA BCT CE contains chemical additives, it is important to avoid possible backflow from the tube.
 To guard against backflow, observe the following precautions:
 - Keep patient's arm in the downward position during the collection procedure.
 - Hold the tube with the stopper in the uppermost position so that the tube contents do not touch the stopper or the end of the needle during sample collection.
 - Release tourniquet once blood starts to flow in the tube, or within 2 minutes of application.
- Follow recommendations for order of draw outlined in CLSI H3-A6⁶.
- Fill tube completely.
- Remove tube from adapter and immediately mix by gentle inversion 8 to 10 times. Inadequate or delayed mixing may result in inaccurate test results. One inversion is a complete turn of the wrist, 180 degrees, and back per the figure below:

- After collection, transport and store tubes within the recommended temperature range.
- Perform extraction in accordance with instrument manufacturer's instructions. For optimal results, please follow the directions for cell-free plasma DNA and cellular genomic DNA extraction.

CELL-FREE PLASMA DNA AND CELLULAR GENOMIC DNA EXTRACTION

- Extraction of cell-free plasma DNA and cellular genomic DNA can be accomplished using most commercially available kits.
- For optimal results, include a Proteinase K treatment step (≥ 30 mAU/ml digest) at 60° C in the presence of chaotropic salts for 1 hour when extracting cell-free DNA and for 2 hours when extracting cellular genomic DNA.

Note:

- Cell-Free DNA BCT CE does not dilute blood samples; therefore, no dilution factor correction is necessary to obtain absolute count values.
- As in the case with most clinical laboratory specimens, hemolysis, icterus and lipemia may affect the results obtained on blood samples preserved with Cell-Free DNA BCT CE.

LIMITATIONS

- Unused tubes to be stored between 18-30° C.
- Samples drawn in other anticoagulants or preservatives may cause coagulation in Cell-Free DNA BCT CE.

REFERENCES

- Das, K., Dumais, J., Basiaga, S., and Krzyzanowski, G., (2012). Carbon-13 Nuclear Magnetic Resonance Analysis of Formaldehyde Free Preservatives. *Acta Histochemica*, 115: 481-486.
- Das, K., Fernando, M.R., Basiaga, S., Wigginton, S., and Williams, T. (2013). Effects of a Novel Cell Stabilizing Reagent on DNA Amplification by PCR as Compared to Traditional Stabilizing Reagents. *Acta Histochemica*, 116: 55-60.
- Norton, S.E., Lechner, J.M., Williams, T., and Fernando, M.R., (2013). A Stabilizing Reagent Prevents Cell-free DNA Contamination by Cellular DNA in Plasma during Blood Sample Storage and Shipping as Determined by Digital PCR. *Clinical Biochemistry*, 46: 1561-1565.
- Norton, S.E., Luna, K.K., Lechner, J.M., Qin, J., and Fernando, M.R., (2013). A New Blood Collection Device Minimizes Cellular DNA Release During Sample Storage & Shipping When Compared to a Standard Device." *Journal of Clinical Laboratory Analysis*, 27: 305-311.
- Qin, J., Alt, J.R., Hunsley, B., Williams, T., and Fernando, M.R. (2013). Stabilization of Circulating Tumor Cells in Blood Using a Collection Device with a Preservative Reagent. *Cancer Cell International*, 14:23.
- Clinical and Laboratory Standards Institute. H3-A6, *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Sixth Edition*.

ORDERING INFORMATION

Please call our Customer Service Department at 402-333-1982 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF HARMONIZED SYMBOLS

EC REP Authorized Representative in the European Community	LOT Batch Code	Biological Risk	REF Catalog Number	Use By
IVD In Vitro Diagnostic Medical Device	Manufacturer	Consult Instructions For Use	Temperature Limitation	Do Not Re-use
Glossary of symbols may contain symbols not used in the labeling of this product.				

See www.streck.com/patents for patents that may be applicable to this product.

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BIBLIOGRAPHIE

1. Das, K., Dumais, J., Basiaga, S., et Krzyzanowski, G., (2012). Carbon-13 Nuclear Magnetic Resonance Analysis of Formaldehyde Free Preservatives. *Acta Histochemica*, 115: 481-486.
2. Das, K., Fernando, M.R., Basiaga, S., Wigginton, S., et Williams, T. (2013). Effects of a Novel Cell Stabilizing Reagent on DNA Amplification by PCR as Compared to Traditional Stabilizing Reagents. *Acta Histochemica*, 116: 55-60.
3. Norton, S.E., Lechner, J.M., Williams, T., et Fernando, M.R., (2013). A Stabilizing Reagent Prevents Cell-free DNA Contamination by Cellular DNA in Plasma during Blood Sample Storage and Shipping as Determined by Digital PCR. *Clinical Biochemistry*, 46: 1561-1565.
4. Norton, S.E., Luna, K.K., Lechner, J.M., Qin, J., et Fernando, M.R., (2013). A New Blood Collection Device Minimizes Cellular DNA Release During Sample Storage & Shipping When Compared to a Standard Device." *Journal of Clinical Laboratory Analysis*, 27: 305-311.
5. Qin, J., Alt, J.R., Hunsley, B., Williams, T., et Fernando, M.R. (2013). Stabilization of Circulating Tumor Cells in Blood Using a Collection Device with a Preservative Reagent. *Cancer Cell International*, 14:23.
6. Clinical and Laboratory Standards Institute. H3-A6, *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Sixième édition*.

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QUELLENANGABEN

1. Das, K., Dumais, J., Basiaga, S. und Krzyzanowski, G., (2012). Carbon-13 Nuclear Magnetic Resonance Analysis of Formaldehyde Free Preservatives. *Acta Histochemica*, 115: 481-486.
2. Das, K., Fernando, M.R., Basiaga, S., Wigginton, S. und Williams, T. (2013). Effects of a Novel Cell Stabilizing Reagent on DNA Amplification by PCR as Compared to Traditional Stabilizing Reagents. *Acta Histochemica*, 116: 55-60.
3. Norton, S.E., Lechner, J.M., Williams, T. und Fernando, M.R., (2013). A Stabilizing Reagent Prevents Cell-free DNA Contamination by Cellular DNA in Plasma during Blood Sample Storage and Shipping as Determined by Digital PCR. *Clinical Biochemistry*, 46: 1561-1565.
4. Norton, S.E., Luna, K.K., Lechner, J.M., Qin, J. und Fernando, M.R., (2013). A New Blood Collection Device Minimizes Cellular DNA Release During Sample Storage & Shipping When Compared to a Standard Device." *Journal of Clinical Laboratory Analysis*, 27: 305-311.
5. Qin, J., Alt, J.R., Hunsley, B., Williams, T. und Fernando, M.R. (2013). Stabilization of Circulating Tumor Cells in Blood Using a Collection Device with a Preservative Reagent. *Cancer Cell International*, 14:23.
6. Clinical and Laboratory Standards Institute. H3-A6, *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Sixth Edition*.

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BIBLIOGRAFIA

1. Das, K., Dumais, J., Basiaga, S., and Krzyzanowski, G., (2012). Carbon-13 Nuclear Magnetic Resonance Analysis of Formaldehyde Free Preservatives. *Acta Histochemica*, 115: 481-486.
2. Das, K., Fernando, M.R., Basiaga, S., Wigginton, S., and Williams, T. (2013). Effects of a Novel Cell Stabilizing Reagent on DNA Amplification by PCR as Compared to Traditional Stabilizing Reagents. *Acta Histochemica*, 116: 55-60.
3. Norton, S.E., Lechner, J.M., Williams, T., and Fernando, M.R., (2013). A Stabilizing Reagent Prevents Cell-free DNA Contamination by Cellular DNA in Plasma during Blood Sample Storage and Shipping as Determined by Digital PCR. *Clinical Biochemistry*, 46: 1561-1565.
4. Norton, S.E., Luna, K.K., Lechner, J.M., Qin, J., and Fernando, M.R., (2013). A New Blood Collection Device Minimizes Cellular DNA Release During Sample Storage & Shipping When Compared to a Standard Device." *Journal of Clinical Laboratory Analysis*, 27: 305-311.
5. Qin, J., Alt, J.R., Hunsley, B., Williams, T., and Fernando, M.R. (2013). Stabilization of Circulating Tumor Cells in Blood Using a Collection Device with a Preservative Reagent. *Cancer Cell International*, 14:23.
6. Clinical and Laboratory Standards Institute. H3-A6, *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Sixth Edition*.

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